

NEW

Info 23.04

"Power Suction" dust removal



Main Switch
Vary infinitely

Plugsocket for castsaw



Hepa Filter

Cast saw with suction for professional use

Bäramed

Instrumente GmbH
Enzianstraße 1
72477 Schwenningen
Tel.: 07579/1730 Fax: 2165
www.Baeramed.de
baeramed@t-online.de



Cast saw "Elektronic Power"



Art.-Nr. V-40100
230V~, 50Hz
Power 500Watt
infinite variable 6500 - 24000 oscillations
5 meter EURO-plug
Plaster saw complete.
1 saw blade Ø50mm for white plaster
1 saw blade Ø65mm for white plaster
1 pair wrenches for saw blade changing
1 instruction manual

This versatile and compact vacuum cleaner with the HEPA filter technology is developed especially for extraction of sawing dust and granulate from ceramics, plastics, metals and hard plasters. The enormous suction power of the turbine is 1450 Watt by an air flow rate of 62 l/s. The vacuum cleaner can be regulated between 500 and 1450 Watt. From the exhaust 99,97% of all particles which are larger than 0,3 micron are filtered. For rational use, the vacuum cleaner is designed with a special outlet for electrical equipment. In combination with the Cast saw Elektronic Power the standby mode automatically takes over the switching on and off of the vacuum function. The cast saw can be placed on the tray. The vacuum cleaner is a robust construction with a special 6 wheel device.
Technical specifications:
230/110 Volt, 50/60 Hz, 1450 Watt



V-308050 saw blade 50mm for hard plaster (blade matt)
V-306050 saw blade 50mm for white plaster (blade gloss)



V-308065 saw blade 65mm for hard plaster (blade matt)
V-306065 saw blade 65mm for white plaster (blade gloss)