

# Neu Bäramed Instrumente GmbH

Enzianstr.1  
72477 Schwenningen  
Tel: +(0) 7579 / 1730 Fax: 2165  
[www.baeramed.de](http://www.baeramed.de)

## Cast Saw „Elektronic Power“



Technical details / specifications  
see overleaf

<b>Plaster Saw „Elektronic Power“</b>	
<b>Art.-No. V-401000</b>	<b>Art.-No. V-501000</b>
<b>230 V~, 50 Hz</b>	<b>110 V~, 60 Hz</b>
<b>Power 500 Watt</b>	<b>Power 500 Watt</b>
<b>6500 – 24000 oscillations, infinitely variable</b>	<b>6500 – 24000 oscillations, infinitely variable</b>
<b>5 meter cable with EURO-plug</b>	<b>3 meter cable wiht USA-plug</b>
Plaster saw complete with carton 1 piece sawblade Ø50 mm for normal plaster 1 piece sawblade Ø65 mm for normal plaster 1 pair wrenches for sawblade changing 1 piece instruction manual	

## Cast Saw „Elektronic Power“

- ✓ **500 Watt power**
- ✓ **6500 bis 24000 oscillations, infinitely variable**
- ✓ **Fullautomatical power-compensation above the whole oscillations field**
- ✓ **Slim motor edition with silent motor 52 – 77 dB(A) in relation to the selected oscillations**
- ✓ **Latest electronic system with various protection and security functions**
- ✓ **Special positioning of the motor, protected by plaster dust**
- ✓ **Robust construction for highest loading capacity by continuous operation**
- ✓ **Safety ON / OFF switch**
- ✓ **Weight only 1,25 kg without cable (1,5 kg with cable)**
- ✓ **Vibrations- and maintenance-free**
- ✓ **Regulation from outside of the manual operating area**
- ✓ **5 meter cable with EU-connection, 230 Volt version or 3 meter cable with US-connection, 110 Volt version**
- ✓ **Double protection isolation according to class II, medical use of typ BF**
- ✓ **Radio and TV suppressed**
- ✓ **CE conformity**
- ✓ **EN 60601-1, EN 60601-1-2, IEC 60601-1, IEC 60601-1-2 inspected**